

Einwegprodukte / Single use



Einweg Doppel-Reinigungsbürste mit Kugelkopf zur Kanalschonung

Disposable double-ended cleaning brush with ball head for channel protection

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Bürste Ø brush	unsteril non-sterile	1 VE pack. unit
FE-D1810WA	1,5 x 1000 mm	2,0 mm	3 mm	⊗	50 Stk./pcs.
FE-D1822WA	1,5 x 2200 mm	2,0 mm	3 mm	⊗	50 Stk./pcs.
FE-D1826WA Enteroskopie	1,5 x 2600 mm	2,0 mm	3 mm	⊗	50 Stk./pcs.
FE-D2410WA	2,4 x 1000 mm	2,8 - 4,2 mm	5 mm	⊗	50 Stk./pcs.
FE-D2422WA*	2,4 x 2200 mm	2,8 - 4,2 mm	5 mm	⊗	50 Stk./pcs.

* auch ohne Kugelköpfe erhältlich / also without ball heads available



Einweg Doppel-Reinigungsbürsten-Set beinhaltet Doppelreinigungsbürste Ø 5 x 2200 mm mit Kugelkopf und Doppel-Ventilreinigungsbürste Ø 11,5 / 5 x 150 mm

Disposable double-ended cleaning brush set contents double-ended cleaning brush Ø 5 x 2200 mm with ball head and double-ended valve cleaning brush Ø 11,5 / 5 x 150 mm

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Bürsten Ø brushes	unsteril non-sterile	1 VE pack. unit
FE-D2422WACR	2,4 x 2200 mm 150 mm	2,8 mm	5 mm 11,5 mm	⊗	100 Sets