


Polypektomie / Polypectomy



Einweg Polypektomieschlinge, asymmetrisch, multifiler Draht, Olympus HF Anschluss

Disposable polypectomy snare, asymmetrical, multifilament wire, Olympus HF connection

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Schlinge Ø snare	steril sterile	1 VE pack. unit
FE-D2423DB25	2,4 x 2300 mm	2,8 mm	25 mm	Ⓢ	10 Stk./pcs.



Einweg Polypektomieschlinge, oval, multifiler Draht, Olympus HF Anschluss, stufenlos verstellbar

Disposable polypectomy snare, oval, multifilament wire, Olympus HF connection, infinitely variable

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Schlinge Ø snare	steril sterile	1 VE pack. unit
FE-D1818DT10	1,8 x 1800 mm	2,0 mm	10 mm	Ⓢ	10 Stk./pcs.
FE-D1818DT20	1,8 x 1800 mm	2,0 mm	20 mm	Ⓢ	10 Stk./pcs.
FE-D2423DT10	2,4 x 2300 mm	2,8 mm	10 mm	Ⓢ	10 Stk./pcs.
FE-D2423DT20	2,4 x 2300 mm	2,8 mm	20 mm	Ⓢ	10 Stk./pcs.
FE-D2423DT30	2,4 x 2300 mm	2,8 mm	30 mm	Ⓢ	10 Stk./pcs.