



Einweg Fangkörbchen zur Steinextraktion, hexagonal, 4-fach mit Anschlagkopf, Notfall-Lithotripsie geeignet
 Disposable retrieval basket for stone extraction, hexagonal with stop head, suitable for emergency lithotripsy

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Korb Ø basket	steril sterile	1 VE pack. unit
FE-D2320LT020	2,3 mm x 2000 mm	ab 2,8 mm	20 mm	Ⓢ	1 Stk./pcs.
FE-D2320LT025	2,3 mm x 2000 mm	ab 2,8 mm	25 mm	Ⓢ	1 Stk./pcs.
FE-D2320LT030	2,3 mm x 2000 mm	ab 2,8 mm	30 mm	Ⓢ	1 Stk./pcs.
FE-D2320LT035	2,3 mm x 2000 mm	ab 2,8 mm	35 mm	Ⓢ	1 Stk./pcs.



Einweg Lithotriptor-Fangkörbchen zur Steinextraktion, hexagonal, 4-fach mit einziehbarer, atraumatischer Spitze, inkl. Spanngriff

Disposable lithotripsy retrieval basket for stone extraction, hexagonal with retractable atraumatic tip, incl. clamp handle

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Korb Ø basket	steril sterile	1 VE pack. unit
FE-D2340LT020	2,3 mm x 4000 mm	ab 2,8 mm	20 mm	Ⓢ	1 Stk./pcs.
FE-D2340LT025	2,3 mm x 4000 mm	ab 2,8 mm	25 mm	Ⓢ	1 Stk./pcs.
FE-D2340LT030	2,3 mm x 4000 mm	ab 2,8 mm	30 mm	Ⓢ	1 Stk./pcs.
FE-D2340LT035	2,3 mm x 4000 mm	ab 2,8 mm	35 mm	Ⓢ	1 Stk./pcs.