



Biopsie Zangen / Biopsy Forceps



Biopsie Zange, oval, distal verjüngte Spirale mit Metallhandgriff, autoklavierbar 134° C

Biopsy forceps, oval, distally tapered spiral with metal handle, autoclavable 134° C

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	unsteril non-sterile	1 VE pack. unit
FE-2417BT	2,4 x 1700 mm	2,8 mm	o	1 Stk./pcs.
FE-2423BT	2,4 x 2300 mm	2,8 mm	o	1 Stk./pcs.
FE-3023BT	3,0 x 2300 mm	3,4 mm	o	1 Stk./pcs.



Biopsie Zange, oval, mit Dorn, distal verjüngte Spirale mit Metallhandgriff, autoklavierbar 134° C

Biopsy forceps, oval, with needle, distally tapered spiral with metal handle, autoclavable 134° C

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	unsteril non-sterile	1 VE pack. unit
FE-2417BW	2,4 x 1700 mm	2,8 mm	o	1 Stk./pcs.
FE-2423BW	2,4 x 2300 mm	2,8 mm	o	1 Stk./pcs.
FE-3023BW	3,0 x 2300 mm	3,4 mm	o	1 Stk./pcs.