

## Reinigung / Cleaning



### Einweg Doppel-Reinigungsbürste mit Kugelkopf zur Kanalschonung

Disposable double-ended cleaning brush with ball head for channel protection

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Bürste Ø brush	unsteril non-sterile	1 VE pack. unit
FE-D1810WA	1,5 x 1000 mm	2,0 mm	3 mm	⊗	50 Stk./pcs.
FE-D1822WA	1,5 x 2200 mm	2,0 mm	3 mm	⊗	50 Stk./pcs.
FE-D1826WA Enteroskopie	1,5 x 2600 mm	2,0 mm	3 mm	⊗	50 Stk./pcs.
FE-D2410WA	2,4 x 1000 mm	2,8 - 4,2 mm	6 mm	⊗	50 Stk./pcs.
FE-D2422WA	2,4 x 2200 mm	2,8 - 4,2 mm	6 mm	⊗	50 Stk./pcs.



### Einweg Dreifach-Reinigungsbürste mit Kugelkopf zur Kanalschonung

Disposable triple-ended cleaning brush with ball head for channel protection

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Bürste Ø brushe	unsteril non-sterile	1 VE pack. unit
FE-D2422WA3	2,4 x 2200 mm	2,8 - 4,2 mm	6 mm	⊗	50 Stk./pcs.