



Einweg Doppel-Reinigungsbürsten-Set beinhaltet Doppelreinigungsbürste Ø 6 x 2200 mm mit Kugelkopf und Doppel-Ventilreinigungsbürste Ø 11,5 / 6 x 150 mm

Disposable double-ended cleaning brush set contents double-ended cleaning brush Ø 6 x 2200 mm with ball head and double-ended valve cleaning brush Ø 11,5 / 6 x 150 mm

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Bürste Ø brushe	unsteril non-sterile	1 VE pack.unit
FE-D2422WACR	2,4 x 2200 mm 150 mm	2,8 mm	6 mm 11,5 / 6 mm	Ⓜ	50 Sets



Einweg Dreifach-Reinigungsbürsten-Set beinhaltet Dreifachreinigungsbürste Ø 6 x 2200 mm mit Kugelkopf und Doppel-Ventilreinigungsbürste Ø 11,5 / 6 x 150 mm

Disposable triple-ended cleaning brush set contents triple-ended cleaning brush Ø 6 x 2200 mm with ball head and double-ended valve cleaning brush Ø 11,5 / 6 x 150 mm

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Bürste Ø brushe	unsteril non-sterile	1 VE pack.unit
FE-D2422WACR3	2,4 x 2200 mm 150 mm	2,8 mm	3 x 6 mm 11,5 / 6 mm	Ⓜ	50 Sets