


Biopsie Zangen / Biopsy Forceps



Einweg Biopsie Zange gezahnt, oval, mit Dorn, beschichtet

Disposable biopsy forceps serrated, oval, with needle, coated

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	steril sterile	1 VE pack. unit
FE-D2418BZC	2,4 x 1800 mm	2,8 mm	Ⓢ	10 Stk./pcs.
FE-D2423BZC	2,4 x 2300 mm	2,8 mm	Ⓢ	10 Stk./pcs.
FE-D2823BZC Jumbo	2,8 x 2300 mm	3,2 mm	Ⓢ	10 Stk./pcs.



Einweg Hot-Biopsie Zange, oval, Olympus HF Anschluss, beschichtet

Disposable Hot biopsy forceps, oval, Olympus HF connection, coated

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	steril sterile	1 VE pack. unit
FE-D2423DG	2,4 x 2300 mm	2,8 mm	Ⓢ	10 Stk./pcs.