




Einweg Polypektomieschlinge, hexagonal, multifiler Draht, Olympus HF Anschluss

Disposable polypectomy snare, hexagonal, multifilament wire, Olympus HF connection

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Schlinge Ø snare	steril sterile	1 VE pack.unit
FE-D2423DL30	2,4 x 2300 mm	2,8 mm	30 mm	Ⓢ	10 Stk./pcs.



Einweg Polypektomieschlinge, oval, monofiler Draht, Olympus HF Anschluss, stufenlos verstellbar

Disposable polypectomy snare, oval, monofilament wire, Olympus HF connection, infinitely variable

Artikelnummer item no.	Ø und Länge Ø and length	Arbeitskanal working channel	Ø Schlinge Ø snare	steril sterile	1 VE pack.unit
FE-D2423DM10	2,4 x 2300 mm	2,8 mm	10 mm	Ⓢ	10 Stk./pcs.
FE-D2423DM20	2,4 x 2300 mm	2,8 mm	20 mm	Ⓢ	10 Stk./pcs.
FE-D2423DM30	2,4 x 2300 mm	2,8 mm	30 mm	Ⓢ	10 Stk./pcs.