




Einweg Doppel-Reinigungsbürsten-Set beinhaltet Doppelreinigungsbürste Ø 6 x 2200 mm mit Kugelkopf und Doppel-Ventilreinigungsbürste Ø 11,5 / 6 x 150 mm

Disposable double-ended cleaning brush set contents double-ended cleaning brush Ø 6 x 2200 mm with ball head and double-ended valve cleaning brush Ø 11,5 / 6 x 150 mm

| Artikelnummer item no. | Ø und Länge Ø and length | Arbeitskanal working channel | Ø Bürste Ø brushe | unsteril non-sterile | 1 VE pack.unit |
|---------------------------|-----------------------------|---------------------------------|----------------------|-------------------------|-------------------|
| FE-D2422WACR | 2,4 x 2200 mm 150 mm | 2,8 mm | 6 mm 11,5 / 6 mm | ⊗ | 50 Sets |



Einweg Dreifach-Reinigungsbürsten-Set beinhaltet Dreifachreinigungsbürste Ø 6 x 2200 mm mit Kugelkopf und Doppel-Ventilreinigungsbürste Ø 11,5 / 6 x 150 mm

Disposable triple-ended cleaning brush set contents triple-ended cleaning brush Ø 6 x 2200 mm with ball head and double-ended valve cleaning brush Ø 11,5 / 6 x 150 mm

| Artikelnummer item no. | Ø und Länge Ø and length | Arbeitskanal working channel | Ø Bürste Ø brushe | unsteril non-sterile | 1 VE pack.unit |
|---------------------------|-----------------------------|---------------------------------|-------------------------|-------------------------|-------------------|
| FE-D2422WACR3 | 2,4 x 2200 mm 150 mm | 2,8 mm | 3 x 6 mm 11,5 / 6 mm | ⊗ | 50 Sets |